

# Identification of selected cations and anions

# Exercise goal:

- pipetting in practise
- identification of selected cations and anions in solutions

Chemical analysis - a set of activities leading to the determination of the chemical, qualitative and quantitative composition of the test substance.

- Qualitative analysis - which determines the chemical elements and compounds consists of the substance; involves the subsequent separation of ion groups using group reagents and detection of individual ions of a given group after the separation of interfering ions
- quantitative analysis - deals with determining the quantitative composition of the substance, i.e. the content of individual components.

<u>Ca<sup>2+</sup></u>	<u>Mg<sup>2+</sup></u>	<u>NH<sub>4</sub><sup>±</sup></u>	<u>Cu<sup>2+</sup></u>	<u>Fe<sup>3+</sup></u>	<u>Zn<sup>2+</sup></u>	<u>Mn<sup>2+</sup></u>	<u>Co<sup>2+</sup></u>	<u>Na<sup>+</sup></u>
		<u>Cl<sup>-</sup></u>	<u>CO<sub>3</sub><sup>2-</sup></u>	<u>PO<sub>4</sub><sup>3-</sup></u>	<u>SO<sub>4</sub><sup>2-</sup></u>	<u>NO<sub>3</sub><sup>-</sup></u>		

## Cation analysis

Cations are divided into 5 groups. This division is based on the reactions that cations of a given group give with group reagents.

Group I:  $\text{Ag}^+$  (silver (I)),  $\text{Hg}_2^{2+}$  (mercury (I)),  $\text{Pb}^{2+}$  (lead (II)).

**Group reagent: HCl (hydrochloric acid) concentration 2 mol /  $\text{dm}^3$ .**

Group II: **Group reagent :  $\text{H}_2\text{S}$  (hydrogen sulfide) in acid medium (HCl concentration 2 mol/ $\text{dm}^3$  ).** Due to the chemical nature of the respective sulphides, group II cations are divided into two subgroups:

Group IIA:  $\text{Cu}^{2+}$  (copper(II)),  $\text{Hg}^{2+}$  (mercury(II)),  $\text{Pb}^{2+}$  (lead(II)) in diluted solutions,  $\text{Cd}^{2+}$  (cadmium(II)) i  $\text{Bi}^{3+}$  (bismuth(III)).  
**Sulphides of these cations are insoluble in  $(\text{NH}_4)_2\text{S}_2$  i KOH.**  
(ammonium sulfide and potassium hydroxide)

Group IIB:  $\text{Sn}^{2+}$  (tin (II)),  $\text{Sn}^{4+}$  (tin(IV)),  $\text{Sb}^{3+}$  (antimony(III)),  $\text{Sb}^{5+}$  (antimony(V)),  $\text{As}^{3+}$  (arsenic(III)) i  $\text{As}^{5+}$  (arsenic(V)).  
**Sulphides of these cations are soluble in  $(\text{NH}_4)_2\text{S}_2$  i KOH.** (ammonium sulfide and potassium hydroxide)

Group III: Group III cations were divided into two subgroups.

Group IIIA:  $\text{Cr}^{3+}$  (chrome(III)),  $\text{Al}^{3+}$  (aluminium),  $\text{Fe}^{2+}$  (iron(II)) i  $\text{Fe}^{3+}$  (iron(III)).

**Group reagent :**  $\text{NH}_3(\text{aq})$  (amonia solution) in the presence of  $\text{NH}_4\text{Cl}$  (ammonium chloride) (weakly alkaline ammonium buffer).

Group IIIB:  $\text{Mn}^{2+}$  (manganese(II)),  $\text{Zn}^{2+}$  (zinc),  $\text{Co}^{2+}$  (cobalt(II)) i  $\text{Ni}^{2+}$  (nickel(II)). **Group reagent :**  $(\text{NH}_4)_2\text{S}$  (ammonium sulfide )in the presence of ammonium buffer.

Group IV:  $\text{Ca}^{2+}$  (calcium),  $\text{Sr}^{2+}$  (strontium) i  $\text{Ba}^{2+}$  (barium).

**Group reagent** :  $(\text{NH}_4)_2\text{CO}_3$  (ammonium carbonate )in the presence of ammonium buffer.

Group V :  $\text{Mg}^{2+}$  (magnesium),  $\text{K}^+$  (potassium),  $\text{NH}_4^+$  (ammonium ion) i  $\text{Na}^+$  (sodium). **This group has no group reagent.**

# Ion detection $\text{Na}^+$

Measure  $1 \text{ cm}^3$   $0,5 \text{ mmol/dm}^3$  of sodium ion into the tube. In this solution, immerse the previously roasted in the flame of the burner platinum wire embedded in a glass rod. Wetted with sodium ion solution, put wire into flame. The flame turns an intense yellow color.

The intensity of the flame color depends on the ion concentration



Lithium




Potassium




Calcium




Nie zabarwiony płomień



Barium




Strontium



Sodium







### Tested Sample

HCl

white precipitate  
Group I:  $\text{Ag}^+$ ,  $\text{Hg}_2^{2+}$ ,  $\text{Pb}^{2+}$

lack precipitate Group II-V

$\text{H}_2\text{S}$  (HCl)

precipitate black, brown or yellow  
Group IIA:  $\text{Cu}^{2+}$ ,  $\text{Hg}^{2+}$ ,  $\text{Pb}^{2+}$ ,  $\text{Cd}^{2+}$ ,  $\text{Bi}^{3+}$

Group IIB:  $\text{Sn}^{2+}$ ,  $\text{Sn}^{4+}$ ,  $\text{Sb}^{3+}$ ,  $\text{Sb}^{5+}$ ,  $\text{As}^{3+}$ ,  $\text{As}^{5+}$

lack precipitate (group III-V)

### Tested Sample

precipitate brown, green or white  
Group IIIA:  $\text{Cr}^{3+}$ ,  $\text{Al}^{3+}$ ,  $\text{Fe}^{2+}$ ,  $\text{Fe}^{3+}$

lack precipitate group IIIB-V

$(\text{NH}_4)_2\text{S}$  ammonium buffer  $\text{H}_3(\text{aq}) + \text{NH}_4\text{Cl}$

precipitate pink white or black  
group IIIB:  $\text{Mn}^{2+}$ ,  $\text{Zn}^{2+}$ ,  $\text{Co}^{2+}$ ,  $\text{Ni}^{2+}$

lack precipitate group IV-V

### Tested Sample

$(\text{NH}_4)_2\text{CO}_3$  ammonium buffer  $\text{H}_3(\text{aq}) + \text{NH}_4\text{Cl}$

white precipitate  
group IV:  $\text{Ca}^{2+}$ ,  $\text{Sr}^{2+}$ ,  $\text{Ba}^{2+}$

lack precipitate  
group V:  $\text{NH}_4^+$ ,  $\text{K}^+$ ,  $\text{Na}^+$

# Anion analysis

In order to facilitate the analysis of anions and standardize analytical procedures, the most common anions were divided into seven analytical groups from I - VII. The division is based on the precipitation of barium and silver salts of individual anions and testing their solubility in water and nitric acid (V)

Group	Solubility, color of precipitate	Detected Ions
<b>I</b>	Ag <sup>+</sup> white or yellow, insoluble in 2M HNO <sub>3</sub> Ba <sup>2+</sup> - easily soluble in water, no precipitate	Cl <sup>-</sup> , Br <sup>-</sup> , F <sup>-</sup> , CN <sup>-</sup> , SCN <sup>-</sup> , ClO <sup>-</sup> , [Fe(CN) <sub>6</sub> ] <sup>4-</sup> , [Fe(CN) <sub>6</sub> ] <sup>3-</sup>
<b>II</b>	Ag <sup>+</sup> - white, soluble in 2M HNO <sub>3</sub> , slightly soluble in water Ba <sup>2+</sup> - precipitate	S <sup>2-</sup> , CH <sub>3</sub> COO <sup>-</sup> , NO <sub>2</sub> <sup>-</sup> , HCOO <sup>-</sup>
<b>III</b>	Ag <sup>+</sup> - colorless, soluble in HNO <sub>3</sub> , insoluble in water Ba <sup>2+</sup> - white soluble soluble in HNO <sub>3</sub> , soluble in water	SO <sub>3</sub> <sup>2-</sup> , CO <sub>3</sub> <sup>2-</sup> , C <sub>2</sub> O <sub>4</sub> <sup>2-</sup> , HPO <sub>4</sub> <sup>2-</sup>
<b>IV</b>	Ag <sup>+</sup> - yellow, light yellow soluble in HNO <sub>3</sub> , insoluble in water Ba <sup>2+</sup> - white precipitates soluble in HNO <sub>3</sub> , hardly soluble in water	S <sub>2</sub> O <sub>4</sub> <sup>2-</sup> , CrO <sub>4</sub> <sup>2-</sup> , Cr <sub>2</sub> O <sub>7</sub> <sup>2-</sup> , PO <sub>4</sub> <sup>3-</sup> , AsO <sub>3</sub> <sup>3-</sup>
<b>V</b>	Ag <sup>+</sup> and barium salts do not precipitate	NO <sub>3</sub> <sup>-</sup> , ClO <sub>3</sub> <sup>-</sup> , ClO <sub>4</sub> <sup>-</sup> , MnO <sub>4</sub> <sup>-</sup>
<b>VI</b>	Ag <sup>+</sup> - do not precipitate Ba <sup>2+</sup> - white precipitates, insoluble in water	SO <sub>4</sub> <sup>2-</sup> , F <sup>-</sup>
<b>VII</b>	Ag <sup>+</sup> - precipitate yellow precipitates, soluble in HNO <sub>3</sub> Ba <sup>2+</sup> - white precipitates, soluble in HNO <sub>3</sub> . Gelatinous precipitates of silicic or tungsten acids are formed	SiO <sub>3</sub> <sup>2-</sup> , WO <sub>4</sub> <sup>2-</sup>

In the following it belongs perform characteristic reactions that will uniquely identify the anion within the group.

Taking into account the oxidation-reducing properties of the above-mentioned anions, we divide them into four groups.

- Anions with reducing properties. These include : chloride anion ( $\text{Cl}^-$ ), iodine anion ( $\text{I}^-$ ), sulfide anion ( $\text{S}^{2-}$ ), sulfate anion ( $\text{SO}_3^{2-}$ ), oxalate anion ( $\text{C}_2\text{O}_4^{2-}$ ) i thiosulphate anion ( $\text{S}_2\text{O}_3^{2-}$ ).
- Anions with oxidizing properties. The anion with strong oxidizing properties is nitrate anion ( $\text{NO}_3^-$ ).
- Anions with both oxidizing and reducing properties. Nitrite anion belongs to this group ( $\text{NO}_2^-$ ).
- Anions that do not show reducing or oxidizing properties in dilute aqueous solutions. We include anions: carbonate anion ( $\text{CO}_3^{2-}$ ), phosphate anion ( $\text{PO}_4^{3-}$ ), sulfate anion ( $\text{SO}_4^{2-}$ ), silicate anion ( $\text{SiO}_3^{2-}$ ).

# Solubility Table

	NH <sub>4</sub> <sup>+</sup>	Na <sup>+</sup>	K <sup>+</sup>	Mg <sup>2+</sup>	Ca <sup>2+</sup>	Ba <sup>2+</sup>	Gf <sup>3+</sup>	Mn <sup>2+</sup>	Fe <sup>2+</sup>	Fe <sup>3+</sup>	Co <sup>2+</sup>	Ni <sup>2+</sup>	Cu <sup>2+</sup>	Ag <sup>+</sup>	Au <sup>3+</sup>	Zn <sup>2+</sup>	Cd <sup>2+</sup>	Hg <sup>2+</sup>	Al <sup>3+</sup>	Sn <sup>2+</sup>	Pb <sup>2+</sup>	Bi <sup>3+</sup>
OH <sup>-</sup>	R	R	R	NR	TR	R	NR	NR	NR	NR	NR	NR	NR	&&	NR	NR	NR	NR	NR	NR	NR	NR
F <sup>-</sup>	R	R	R	NR	NR	NR	NR	TR	TR	NR	R	TR	NR	R	R	NR	TR	R	TR	R	NR	R
Cl <sup>-</sup>	R	R	R	R	R	R	R	R	R	R	R	R	R	NR	R	R	R	R	R	R	TR	R
Br <sup>-</sup>	R	R	R	R	R	R	R	R	R	R	R	R	R	NR	R	R	R	TR	R	R	TR	R
I <sup>-</sup>	R	R	R	R	R	R	R	R	R	&&	R	R	&&	NR	TR	R	R	NR	R	TR	NR	NR
S <sup>2-</sup>	R	R	R	&&	TR	R	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	&&	NR	NR	NR
SO <sub>3</sub> <sup>2-</sup>	R	R	R	R	NR	NR	R	NR	NR	&&	NR	NR	&&	TR	NR	TR	TR	&&	&&	TR	NR	TR
SO <sub>4</sub> <sup>2-</sup>	R	R	R	R	TR	NR	R	R	R	R	R	R	R	TR	R	R	R	R	R	R	NR	R
NO <sub>2</sub> <sup>-</sup>	R	R	R	R	R	R	R	R	R	NR	R	R	R	TR	&&	R	R	R	R	TR	R	NR
NO <sub>3</sub> <sup>-</sup>	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R	R
PO <sub>4</sub> <sup>2-</sup>	R	R	R	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR	NR
CO <sub>3</sub> <sup>2-</sup>	R	R	R	NR	NR	NR	NR	NR	NR	&&	NR	NR	&&	NR	NR	NR	NR	NR	&&	NR	NR	NR
SiO <sub>3</sub> <sup>2-</sup>	R	R	R	NR	NR	NR	&&	NR	NR	NR	NR	NR	&&	&&	&&	NR	NR	&&	NR	&&	NR	&&
MnO <sub>4</sub> <sup>2-</sup>	R	R	R	R	R	R	R	&&	&&	R	R	R	R	R	&&	R	R	&&	R	&&	R	R
CrO <sub>4</sub> <sup>2-</sup>	R	R	R	R	TR	NR	NR	NR	&&	R	NR	NR	R	NR	&&	TR	NR	TR	NR	NR	NR	NR
octan	R	R	R	R	R	R	R	NR	R	R	R	R	R	R	R	R	R	R	R	R	R	R

Kolor tła odzwierciedla charakterystyczną barwę substancji, należy go traktować umownie.

R - substancja dobrze rozpuszczalna (rozpuszczalność powyżej 1 g w 100 g wody).

TR - substancja o niewielkiej rozpuszczalności, strąca się przy odpowiednim stężeniu roztworu (rozpuszczalność 0,1 - 1 g w 100 g wody).

NR - substancja praktycznie nierozpuszczalna, strąca się z rozcieńczonych roztworów (rozpuszczalność poniżej 0,1 g w 100 g wody).

&& - zachodzą skomplikowane reakcje, lub substancja nie została otrzymana.

### Task 1


The aim of the task is to detect  $\text{Ca}^{2+}$  (calcium) ions

#### Execution

##### a. Test with $\text{C}_2\text{O}_4^{2-}$ (oxalate) ions

Add about  $0.5 \text{ cm}^3$  of  $0.2 \text{ mol/dm}^3$  ammonium oxalate ( $\text{NH}_4\text{C}_2\text{O}_4$ ) to  $0.5 \text{ cm}^3$  of  $\text{Ca}^{2+}$  solution.

Reaction:




The resulting  $\text{CaC}_2\text{O}_4$  precipitate is white, crystalline. If you add several drops of ammonia to a test tube, in which the reaction has occurred, the amount of precipitate is greater. Calcium oxalate is soluble in mineral (inorganic) acids, and insoluble in acetic acid.

##### b. Test with $\text{SO}_4^{2-}$ (sulphate) ions

Add several drops of concentrated  $\text{H}_2\text{SO}_4$  (sulphuric acid) to  $0.5 \text{ cm}^3$  of  $\text{Ca}^{2+}$  solution.

Reaction:



Reaction with sulphuric acid results in the precipitation of white  $\text{CaSO}_4$  (copper sulphate).



## Task 2

The aim of the task is to detect  $\text{NH}_4^+$  (ammonium) ions


### Execution

Add about  $0.5 \text{ cm}^3$  of  $2 \text{ mol/dm}^3 \text{ NaOH}$  to  $0.5 \text{ cm}^3$  of  $\text{NH}_4^+$  solution, and heat the contents of the test tube over a burner.

Reaction:



Observe the change of litmus paper colour and the characteristic odour of escaping ammonia at the mouth of the test tube.





### Task 3

The aim of the task is to detect  $\text{Cu}^{2+}$  (copper) ions

#### Execution

Take two test tubes, pour into each  $0.5 \text{ cm}^3$  of  $\text{Cu}^{2+}$  solution, and then add  $0.5 \text{ cm}^3$  of  $2 \text{ mol/dm}^3 \text{ NaOH}$ . After the precipitation of  $\text{Cu}(\text{OH})_2$ , heat one test tube over a burner, and add  $0.5 \text{ cm}^3$  of concentrated  $\text{NH}_3 \cdot \text{H}_2\text{O}$  to the other test tube.

Reactions:





#### Task 4


The aim of the task is to detect  $\text{Fe}^{3+}$  (iron) ions

Execution

a. Test with potassium thiocyanate

Add  $0.5 \text{ cm}^3$  of  $1 \text{ mol/dm}^3$  KSCN (potassium thiocyanate) to  $0.5 \text{ cm}^3$  of  $\text{Fe}^{3+}$  solution.

Reaction:




Observe the change of solution colour.  $\text{Fe}^{2+}$  ions do not give this reaction.


b. Test with sodium hydroxide

Add  $0.5 \text{ cm}^3$  of  $2 \text{ mol/dm}^3$  NaOH to  $0.5 \text{ cm}^3$  of  $\text{Fe}^{3+}$  solution.

Reaction:



Observe the formation of red-brown  $\text{Fe}(\text{OH})_3$  precipitate, which is insoluble in an excess of the reagent.



### Task 5



The aim of the task is to detect  $\text{Ag}^+$  (silver) ions

#### Execution



##### a. Test with $\text{Cl}^-$ (chlorine) ions

Take two test tubes, pour into each  $0.5 \text{ cm}^3$  of  $\text{Ag}^+$  solution, and then several drops of concentrated HCl (hydrochloric acid). One test tube expose to sunlight and after some time observe photochemical reaction occurring in the test tube. Carefully heat the other test tube and observe the precipitate. Cool the test tube, then add drop by drop  $2 \text{ mol/dm}^3 \text{ NH}_3 \cdot \text{H}_2\text{O}$  until the formation of soluble silver diamine chloride ( $[\text{Ag}(\text{NH}_3)_2]\text{Cl}$ ), and next add drop by drop concentrated  $\text{HNO}_3$  until the reappearance of AgCl precipitate.


$\text{Ag}^+$  ions, occurring in neutral or acid solutions, react with  $\text{Cl}^-$  ions forming white, cheese-like precipitate of silver chloride (AgCl), which decomposes when exposed to light and darkens due to the precipitation of colloidal silver.



AgCl precipitate is insoluble in both cold and hot pure water but in the presence of ammonium ions it forms soluble, colourless, complex ion  $[\text{Ag}(\text{NH}_3)_2]^+$ .




The addition of nitric acid causes the reprecipitation of AgCl.




**b. Test with  $\text{CrO}_4^{2-}$  (chromate) ions**

Add  $0.5 \text{ cm}^3$  of  $0.1 \text{ mol/dm}^3 \text{ K}_2\text{CrO}_4$  to  $0.5 \text{ cm}^3$  of  $\text{Ag}^+$  solution. After the precipitation of  $\text{Ag}_2\text{CrO}_4$ , add several drops of concentrated nitric acid to the test tube. Observe the course of reaction.



Potassium chromate ( $\text{K}_2\text{CrO}_4$ ) causes the precipitation of red-brown silver chromate ( $\text{Ag}_2\text{CrO}_4$ ) from the solutions containing  $\text{Ag}^+$  ions. This reaction has found practical application in the quantitative determination of silver in jewellery because the intensity of colour of formed  $\text{Ag}_2\text{CrO}_4$  depends on the amount of silver in the alloy, from which jewellery has been made.



## Task 6


The aim of the task is to detect  $\text{Pb}^{2+}$  (lead) ions

### Execution

#### a. Test with $\text{Cl}^-$ (chlorine) ions

Add several drops of concentrated HCl to  $0.5 \text{ cm}^3$  of  $\text{Pb}^{2+}$  solution.

Reaction:



The precipitate of lead chloride (II) is formed. Add  $2 \text{ cm}^3$  of distilled water to the test tube and heat until the precipitate is dissolved. After cooling, observe the formed crystals.


$\text{Cl}^-$  ions, added to the solutions containing  $\text{Pb}^{2+}$  ions, cause the precipitation of white, crystalline lead chloride (II), which quickly settle on the test tube bottom.

This precipitate dissolves after heating in the test tube, and after cooling it reprecipitates in the form of needle-like crystals.



Test with heating silver and lead chlorides is important because it is used to easily distinguish AgCl from  $\text{PbCl}_2$  (AgCl is insoluble in hot water).

#### b. Test with chromate

Add  $0.5 \text{ cm}^3$  of  $0.1 \text{ mol/dm}^3 \text{ K}_2\text{CrO}_4$  to  $0.5 \text{ cm}^3$  of  $\text{Pb}^{2+}$  solution. After a few minutes, examine the formed precipitate.



Potassium chromates ( $\text{K}_2\text{CrO}_4$  and  $\text{K}_2\text{Cr}_2\text{O}_7$ ), added to the solutions containing  $\text{Pb}^{2+}$  ions, causes the precipitation of yellow lead chromates (II).



### Task 7


The aim of the task is to detect  $\text{Ba}^{2+}$  (barium) ions

#### Execution

##### a. Test with $\text{SO}_4^{2-}$ (sulphate) ions

Add 0.5 cm<sup>3</sup> of 0.1 mol/dm<sup>3</sup>  $\text{H}_2\text{SO}_4$  to 0.5 cm<sup>3</sup> of  $\text{Ba}^{2+}$  solution.

Reaction:





Sulphuric acid (VI) causes the precipitation of white  $\text{BaSO}_4$ , which is insoluble in water.

##### b. Test with $\text{CO}_3^{2-}$ (carbonate) ions



Add 0.5 cm<sup>3</sup> of 0.2 mol/dm<sup>3</sup>  $\text{Na}_2\text{CO}_3$  to 0.5 cm<sup>3</sup> of  $\text{Ba}^{2+}$  solution.

After the precipitation of  $\text{BaCO}_3$ , add 1 cm<sup>3</sup> of concentrated  $\text{CH}_3\text{COOH}$ .

Reactions:



Carbonate ions, added to the solutions containing  $\text{Ba}^{2+}$  ions, cause the precipitation of white barium carbonate (II), which is easily soluble in acetic acid. If the precipitate is not dissolved, you can add 1 cm<sup>3</sup> of  $\text{CH}_3\text{COOH}$  more.



### Task 8


The aim of the task is to detect  $\text{Cl}^-$  (chloride) ions

#### Execution

##### a. Test with $\text{Ag}^+$ (silver) ions

Add drop by drop about  $0.5 \text{ cm}^3$  of  $0.1 \text{ mol/dm}^3 \text{ AgNO}_3$  to  $0.5 \text{ cm}^3$  of  $\text{Cl}^-$  solution. Dissolve the resulting precipitate adding  $2 \text{ mol/dm}^3 \text{ NH}_3 \cdot \text{H}_2\text{O}$  drop by drop. Then add several drops of concentrated  $\text{HNO}_3$  in order to cause the reprecipitation of  $\text{AgCl}$ .

Reactions:





##### b. Test with potassium permanganate

Add  $0.5 \text{ cm}^3$  of  $0.1 \text{ mol/dm}^3 \text{ KMnO}_4$  and several drops of concentrated  $\text{H}_2\text{SO}_4$  to  $0.5 \text{ cm}^3$  of  $\text{Cl}^-$  solution.

Reaction:



Potassium permanganate, as strong oxidant, is reduced in acid environment (the solution becomes discoloured), while  $\text{Cl}^-$ , as weak reducer, escapes in the form of gas with a characteristic odour.






### Task 9

The aim of the task is to detect  $\text{CO}_3^{2-}$  ions

#### Execution


##### a. Test with strong acid

Add drop by drop 2 mol/dm<sup>3</sup>  $\text{H}_2\text{SO}_4$  to 0.5 cm<sup>3</sup> of  $\text{CO}_3^{2-}$  solution. Observe the course of reaction.




##### b. Test with $\text{Ba}^{2+}$ ions

Add 0.5 cm<sup>3</sup> of 0.1 mol/dm<sup>3</sup>  $\text{BaCl}_2$  to 0.5 cm<sup>3</sup> of  $\text{CO}_3^{2-}$  solution. Reaction:



White precipitate of barium carbonate is soluble in acetic, nitric and hydrochloric acids, and insoluble in sulphuric acid.




### Task 10

The aim of the task is to detect  $\text{PO}_4^{3-}$  ions

#### Execution

Add 5 drops of concentrated  $\text{HNO}_3$  and  $0.5 \text{ cm}^3$  of  $(\text{NH}_4)_2\text{MoO}_4$  to  $0.5 \text{ cm}^3$  of  $\text{PO}_4^{3-}$  solution. Heat the test tube over a burner.

Reaction:



### Task 11


The aim of the task is to detect  $\text{SO}_4^{2-}$  ions

#### Execution

##### a. Test with $\text{Pb}^{2+}$ (lead) ions

Add 0.5 cm<sup>3</sup> of 2 mol/dm<sup>3</sup>  $\text{Pb}(\text{CH}_3\text{COO})_2$  to 0.5 cm<sup>3</sup> of  $\text{SO}_4^{2-}$  solution.

Reaction:




The resulting white precipitate of  $\text{PbSO}_4$  is insoluble in inorganic acids, except concentrated sulphuric acid (VI), and it is also insoluble in concentrated NaOH and KOH.

##### b. Test with $\text{Ba}^{2+}$ (barium) ions

Add 0.5 cm<sup>3</sup> of 1 mol/dm<sup>3</sup> barium chloride ( $\text{BaCl}_2$ ) to 0.5 cm<sup>3</sup> of  $\text{SO}_4^{2-}$  solution. Reaction:



The resulting white precipitate of  $\text{BaSO}_4$  is insoluble in HCl and HNO<sub>3</sub> even after heating, while it is soluble in concentrated  $\text{H}_2\text{SO}_4$ .



## Task 12


The aim of the task is to detect  $\text{C}_2\text{O}_4^{2-}$  ions

### Execution

#### a. Test with $\text{CaCl}_2$

Add drop by drop  $0.5 \text{ cm}^3$  of  $1 \text{ mol/dm}^3 \text{ CaCl}_2$  to  $0.5 \text{ cm}^3$  of  $\text{C}_2\text{O}_4^{2-}$  solution.

Reaction:



Calcium chloride causes the precipitation of white crystalline precipitate, which is soluble in mineral acids, and insoluble in acetic acid.



#### b. Test with potassium permanganate

Add  $0.5 \text{ cm}^3$  of  $0.1 \text{ mol/dm}^3 \text{ KMnO}_4$  and  $0.5 \text{ cm}^3$  of concentrated  $\text{H}_2\text{SO}_4$  to  $0.5 \text{ cm}^3$  of  $\text{C}_2\text{O}_4^{2-}$  solution.

Reaction:



Potassium permanganate becomes discoloured in acid solution containing oxalate ions ( $\text{C}_2\text{O}_4^{2-}$ ). Characteristic bubbles of  $\text{CO}_2$  are formed during this reaction.



# Disinfection (decontamination)

Disinfection is a process that is designed to kill actively growing and vegetative microbial microorganisms to prevent infection.

Disinfectants - chemical substances that destroy (kill) all microorganisms and their spore forms, also preventing the reproduction of microorganisms. Aseptic conditions can be achieved with disinfectants.

There are two synonymous terms that should be distinguished: **antiseptics and aseptics.**

# Disinfection (decontamination)

- **Antiseptics** are the destruction of microbes on tissues with antiseptics substances
- **Aseptics** are all proceedings preventing infection with pathogenic microorganisms.

Physical and chemical methods are used in aseptics. Aseptic conditions are an environment devoid of any pathogenic life forms.

# Disinfection (decontamination)


- **Antiseptic preparations** includes substances used to disinfect tissues and prevent infections on the surface of the skin and mucous membranes.
- **Antiseptics** is any chemical compound (or element, e.g. fluorine, ozone, chlorine) that destroys microorganisms and prevents their reproduction.
- **disinfectants** these are substances used to disinfect rooms, objects (instruments, sanitary equipment, furniture, walls, floors), or fragments of the natural environment (e.g. soil, watercourses)




1 <b>H</b> 1.008 Hydrogen																	2 <b>He</b> 4.003 Helium				
3 <b>Li</b> 6.941 Lithium	4 <b>Be</b> 9.012 Beryllium															5 <b>B</b> 10.811 Boron	6 <b>C</b> 12.011 Carbon	7 <b>N</b> 14.007 Nitrogen	8 <b>O</b> 15.999 Oxygen	9 <b>F</b> 18.998 Fluorine	10 <b>Ne</b> 20.180 Neon
11 <b>Na</b> 22.990 Sodium	12 <b>Mg</b> 24.305 Magnesium															13 <b>Al</b> 26.982 Aluminium	14 <b>Si</b> 28.086 Silicon	15 <b>P</b> 30.974 Phosphorus	16 <b>S</b> 32.065 Sulfur	17 <b>Cl</b> 35.453 Chlorine	18 <b>Ar</b> 39.948 Argon
19 <b>K</b> 39.098 Potassium	20 <b>Ca</b> 40.078 Calcium	21 <b>Sc</b> 44.956 Scandium	22 <b>Ti</b> 47.867 Titanium	23 <b>V</b> 50.942 Vanadium	24 <b>Cr</b> 51.996 Chromium	25 <b>Mn</b> 54.938 Manganese	26 <b>Fe</b> 55.845 Iron	27 <b>Co</b> 58.933 Cobalt	28 <b>Ni</b> 58.693 Nickel	29 <b>Cu</b> 63.546 Copper	30 <b>Zn</b> 65.38 Zinc	31 <b>Ga</b> 69.723 Gallium	32 <b>Ge</b> 72.630 Germanium	33 <b>As</b> 74.922 Arsenic	34 <b>Se</b> 78.96 Selenium	35 <b>Br</b> 79.904 Bromine	36 <b>Kr</b> 83.80 Krypton				
37 <b>Rb</b> 85.468 Rubidium	38 <b>Sr</b> 87.62 Strontium	39 <b>Y</b> 88.906 Yttrium	40 <b>Zr</b> 91.224 Zirconium	41 <b>Nb</b> 92.906 Niobium	42 <b>Mo</b> 95.94 Molybdenum	43 <b>Tc</b> 98 Technetium	44 <b>Ru</b> 101.07 Ruthenium	45 <b>Rh</b> 102.91 Rhodium	46 <b>Pd</b> 106.42 Palladium	47 <b>Ag</b> 107.87 Silver	48 <b>Cd</b> 112.41 Cadmium	49 <b>In</b> 114.82 Indium	50 <b>Sn</b> 118.71 Tin	51 <b>Sb</b> 121.76 Antimony	52 <b>Te</b> 127.60 Tellurium	53 <b>I</b> 126.90 Iodine	54 <b>Xe</b> 131.29 Xenon				
55 <b>Cs</b> 132.91 Cesium	56 <b>Ba</b> 137.33 Barium	57 - 71 <b>La - Lu</b>	72 <b>Hf</b> 178.49 Hafnium	73 <b>Ta</b> 180.95 Tantalum	74 <b>W</b> 183.84 Tungsten	75 <b>Re</b> 186.21 Rhenium	76 <b>Os</b> 190.23 Osmium	77 <b>Ir</b> 192.22 Iridium	78 <b>Pt</b> 195.08 Platinum	79 <b>Au</b> 196.97 Gold	80 <b>Hg</b> 200.59 Mercury	81 <b>Tl</b> 204.38 Thallium	82 <b>Pb</b> 207.2 Lead	83 <b>Bi</b> 208.98 Bismuth	84 <b>Po</b> 209 Polonium	85 <b>At</b> 210 Astatine	86 <b>Rn</b> 222 Radon				
87 <b>Fr</b> 223 Francium	88 <b>Ra</b> 226 Radium	89 <b>Ac</b> 227 Actinide	90 <b>Th</b> 232.04 Thorium	91 <b>Pa</b> 231.04 Protactinium	92 <b>U</b> 238.03 Uranium																



Bulk biological elements



Trace elements believed to be essential for bacteria, plants or animals



Possibly essential trace elements for some species

# Calcium (Ca)

- Calcium is an essential component of bone and cartilage
- Calcium is essential for the normal clotting of blood, by stimulating the release of thromboplastin from the blood platelets
- Calcium is an activator for several key enzymes, including pancreatic lipase, acid phosphatase, cholinesterase and succinic dehydrogenase
- Calcium stimulates muscle contraction (ie. promotes muscle tone and normal heart beat)
- Calcium, in conjunction with lipids, plays a key role in the regulation of the permeability of cell membranes
- Calcium occurs in the form of calcium phosphate in bones
- The ionic form of  $\text{Ca}^{2+}$  is present in extracellular fluids

# Manganese

**Potassium permanganate  $\text{KMnO}_4$**  belongs to the group of oxidizing antiseptics. Under the influence of organic compounds, e.g. proteins, it is reduced and releases oxygen that destroys bacteria, fungi and protozoa. It dissolves easily in water. Depending on the concentration, aqueous solutions are pink or dark purple.

- It works astringent, anti-inflammatory, anti-inflammatory, bactericidal, bacteriostatic, fungicidal, fungistatic, antiviral and antiprotozoal.
- Concentrated solutions are coloring and corrosive on the skin. Under the influence of body secretions and enzymes, it gradually turns into manganese dioxide  $\text{MnO}_2$
- Potassium permanganate neutralizes venoms invertebrates and vertebrates.
- In antiseptics are used 0.5-4% solutions. They are used for decontaminating wounds and burns.



# Iodine

**Iodine It is very hard soluble in water, easily soluble in alcohols and acetone.**

Easily soluble in aqueous potassium iodide solution

iodine has antiseptic properties is used to disinfect intact skin. It has bactericidal, fungicidal and virucidal effect, irritates tissues, causes surface necrosis (necrosis) of smeared tissue.

Iodine is used in the form of a 3-10% solution in 70% or 95% ethyl alcohol. Iodine reacts with the amino groups of proteins, inhibiting their activity. It also releases oxygen from aqueous solutions, which complements the sanitizing effect.

- To disinfect the skin - wash with iodine or Lugol's solution. Do not apply undiluted directly to wounds and burns. Undiluted, it is suitable for disinfecting only around wounds, burns and intact skin.





# Iodine

- Iodine is an essential constituent of the thyroid hormones
- Iodine does not accumulate in the body, therefore it must be constantly supplied with the diet
- Iodine deficiency causes hypothyroidism, which has serious metabolic consequences leading to a slowing of vital functions and all cellular changes.



# Factors affecting the effectiveness of biocidal preparations

- contact time of the preparation with microorganisms (+++)
- concentration of the preparation(++)
- temperature (++)
- water hardness and environmental pH(+/-)
- number of microorganisms(--)
- protein load (large amount of organic substance in the environment: food residues, feces)(--)
- the presence of surfactants.(++)

Of the factors outlined above, the duration and concentration of the preparation are undoubtedly the most important