

Buffer solutions

The aim:

- to familiarize with the properties of solutions
- to familiarize with the mechanism of buffers action and the properties of buffer solutions

The mechanism of action of acetate buffer

Acetic acid present in the buffer solution is practically in undissociated form and acts as proton donor. While, sodium acetate has the ion form: Na^+ and CH_3COO^- , wherein CH_3COO^- ions (as a strong base) are proton acceptors. The proton donor - CH_3COOH (the acid according to the Brønsted–Lowry theory) protects the solution from pH changes during adding the base to the system. Added excess of the OH^- ions turns to the undissociated water:



The proton acceptor - CH_3COO^- ions (the base according to the Brønsted–Lowry theory) protects the solution from pH changes during adding small amount of acid to the system (hydrogen ions, oxonium ions - H_3O^+):



During the reaction, poorly dissociated acetic acid is created and pH of the solution practically does not change. In this way, buffer system works on both sides - counteracts the increasing and decreasing pH value.

The mechanism of action of phosphate buffer

Sodium dihydrogen phosphate (V) present in the buffer solution has an ion form – sodium ions (Na^+) and dihydrogen phosphate ions (H_2PO_4^-) and acts as proton donor. While, sodium hydrogen phosphate (V) consists of: sodium ions (Na^+) and hydrogen phosphate ions (HPO_4^{2-}). HPO_4^{2-} ions as a base are proton acceptors. The proton donor: H_2PO_4^- (the acid according to the Brønsted–Lowry theory) protects the solution from pH changes during adding the base to the system. Added excess of the OH^- ions turns to the undissociated water:



The proton acceptor - HPO_4^{2-} ions (the base according to the Brønsted–Lowry theory) protects the solution from pH changes during adding small amount of acid to the system (hydrogen ions, oxonium ions - H_3O^+):



During the reaction dihydrogen phosphate ion is created and pH of the solution practically does not change.

The mechanism of action of an indicator

To precisely conduct the neutralization reaction, it is necessary to use the indicators. Because the indicator changes its color depending on the pH of the solution. Even a very small amount of acid or base causes an immediate color change of the indicator. Change the color of indicator is caused by its dissociation or withdrawal of the dissociation rate.

The conclusion is that the **color of the solution depends on the color of ions or the color of the undissociated indicator molecules**. Therefore, the indicators that we generally use are weak acids or weak organic bases.

The mechanism of action of an indicator

Conventional indicators can be divided into:

- 1) monochromatic
- 2) dichromatic

Each indicator used in neutralization reaction is characterised by so-called the range of color change depending on pH.


Example of dichromatic indicator:

Bromothymol blue changes the color between pH 6.2 and 7.6. It means that this indicator in the solutions of pH = 6.2 or less is **yellow**. If the solution has pH 7.6 or higher, after adding bromothymol blue it becomes **blue**. Whereas, the solutions of pH within the range 6.2 – 7.6, under the influence of bromothymol are **green**. Green is the range of color change of bromothymol blue.




Indicator	pH range for color change
malachite green	0,0 - 2,0
brilliant green	0,0 - 2,6
methyl green	0,1 - 2,3
methyl violet	0,1 - 2,7
cresol red	0,2 - 1,8
thymol blue	1,2 - 2,8
<u>dinitrophenol</u>	2,4 - 4,0
methyl yellow	2,9 - 4,0
methyl orange	3,1 - 4,4
bromocresol green	3,8 - 5,4
methyl red	4,2 - 6,3
litmus	4,5 - 8,3
bromocresol red	5,2 - 6,8
bromothymol blue	6,2 - 7,6
phenol red	6,4 - 8,0
cresol red	7,2 - 8,8
phenolphthalein	8,3 - 10,0
thymolphthalein	9,3 - 10,5
<u>tropeolin O</u>	11,0-13,0

Selected indicators

Bromocresol green







Bromocresol blue



PHENOLPHTHALEIN

structure



pH

< 0

0 – 8.2

8.2 – 12.0

>12.0

color


orange

colorless

from pink
to fuchsia

colorless

PHENOLPHTHALEIN



Solutions

1. True solution
2. Colloidal solution
3. Suspension

True solution

True Solution is a homogeneous mixture of two or more substances in which substance dissolved (solute) in solvent has the particle size of less than 10^{-9} m. Particles of true solution cannot react with each other.

Conventionally, it is assumed that, the ingredient which is in greater amount – is a solvent (dispersion phase) and compound in smaller amount – is a dissolved substance (dispersed phase).

Molecules of the dissolved substance are surrounded by solvent molecules – solvation (if the solvent is water - hydration).

Water solutions are a fundamental component of the living matter and the environment in which there is life.

Colloidal system

Colloidal solution has the molecules so small that it does not form a suspension, but also too large to form a true solution. Colloidal particles are visible under the electron microscope, we can separate them in ultracentrifuges.

Example: blood plasma of animals

Colloids show some unique properties:

- Brownian motion
- Tyndall effect

Suspensions

- These are the heterogenous and biphasic systems, in which the diameter of particles of the substance dispersed in the solvent exceeds 100 nm.
- Examples: soup, mud, water in the lake, etc.

Solubility

Solubility - is the maximum number of grams of a substance which can be dissolved in 100 g of solvent (Temp., Pres. = Const).

Solubility of the substance is equal to the concentration of the saturated solution.

Substance solubility depends on:

- a. type of the dissolved substance (solute)
- b. type of the solvent
- c. temperature (solubility of solids in liquids increases with temperature, and gases decreases)
- d. pressure (solubility of gases in liquids increases with pressure)

Important processes for the living organism

Osmosis – is the spontaneous movement of solvent molecules through a semi-permeable membrane into a region of higher solute concentration.

Dialysis – is the movement of solute particles through a semi-permeable membrane.

Semi-permeable membranes – permeate the molecules of water and other solvents, while they do not permeate many molecules of dissolved substances, especially macromolecules.

It does not mean that naturally occurring membranes are only the mechanical sieves that permeable small molecules and retain larger particles.

They **selectively permeate** one group of molecules, do not allowing other molecules (even of the same size) to cross the membrane, at the same time.

Task 1

The purpose: to prepare and to calculate the percentage and molar concentration of the copper sulfate solutions.

Procedure:

1. Add the following volume of the CuSO_4 solution and distilled water to 5 tubes (I, II, III, IV, V) according to the table below:

Tube No.	CuSO_4 cm^3	H_2O cm^3	Concentration	
			%	mol/dm^3
I	2	1		
II	1	1		
III	1	2		
IV	0.5	4.5		
V	0.1	0.9		

2. Calculate the final concentration (percentage and molar) of the copper sulfate solutions.

Concentration of stock solution was 1%.

Molar mass of $\text{CuSO}_4 = 159.61 \text{ g/mol}$

Task 2

The purpose: to prepare the buffer solutions of a selected pH.

Procedure:

1. Add 0.1 mol/dm³ acetic acid and 0.1 mol/dm³ sodium acetate to 5 tubes (as described in the table below):

Tube No.	CH ₃ COOH 0.1 mol/dm ³ cm ³	CH ₃ COONa 0.1 mol/dm ³ cm ³	Calculated pH value (for 18°C)	pH value measured potentiometri- cally
1	9	1		
2	7	3		
3	5	5		
4	3	7		
5	1	9		

2. Precisely mix the tubes.

3. Calculate the pH of the buffer solutions from each tube using the **Henderson-Hasselbalch equation**.

4. Measure the pH by pH-meter to confirm the accuracy of the calculations.

Task 3

The purpose: to observe the effect of the dilution of buffer solution on its pH.

Procedure:

1. Add 0.1 mol/dm³ acetate buffer (pH=4.7) and distilled water to 5 tubes


Tube No.	Buffer cm ³	H ₂ O cm ³	Dilution	Molarity	Color of the indicator
1	2	0			
2	0.5	1.5			
3	0.2	1.8			
4	0.02	1.98			
5	0	2			

2. Calculate the dilution and molarity of the buffer solutions deriving from each tube.

3. Add 5 drops of bromocresol green to each tube.

4. Consider and answer the question:

What is the color of the indicator at various dilutions of the buffer solution? What does it mean?



Task 4

The purpose: to evaluate the capacity of the acetate buffer.

Procedure:

1. Add 0.1 mol/dm³ sodium acetate, 0.1 mol/dm³ acetic acid and distilled water in appropriate volume to 4 tubes:

Tube No.	CH ₃ COOH 0.1 mol/dm ³ cm ³	CH ₃ COONa 0.1 mol/dm ³ cm ³	H ₂ O cm ³	Dilution	The volume of consumed NaOH (drops)
1	2	2	0		
2	1	1	2		
3	0.5	0.5	3		
4	0	0	4		

2. Calculate the dilution of the buffer solution.

3. Add 5 drops of bromocresol green to each tube.

4. Add 0.1 cm³ of 0.1 mol/dm³ NaOH to tube No. 4.

5. Add by drops 0.1 mol/dm³ NaOH to other tubes, until the color has become the same as in the tube No. 4. **Note the volume of the consumed NaOH in each tube!**

Does the dilution has an influence on the buffer capacity?