

## Body fluids Physiological and pathological urine

Enter obtained results into the table below.

Parameter		Sample	
		Physiological urine	Pathological urine
Color of urine			-
pH of urine			-
Density of urine			-
Chloride ions			-
Phosphate ions			-
Calcium ions			-
Sulphate ions			-
Protein	With sulphosalicylic acid		
	With trichloroacetic acid		
	Reaction of Heller		
Glucose			
Ketone bodies			
Bile pigments	Reaction of Gmelin		
	Reaction of Rosin		

### 1. Physical properties of physiological urine

#### Protocol:

#### 1. Color of urine

Urine of majority of animals is transparent and yellow.

#### 2. pH of urine.

The examination of pH of urine is done by use of paper (strip) indicators. pH of urine of carnivorous is between 5,0-7,0, while herbivorous between 7,4-8,2.

#### 3. Density of urine.

Density of urine is examined with urometer. Pour urine into the cylinder and insert urometer. Read off the scale numerical value which represents density of urine in g/cm<sup>3</sup>. Physiological values are between 1,015-1,060 g/cm<sup>3</sup>

### 2. Chemical properties of physiological urine

#### Protocol.

#### - 1. Detection of Chloride ions

Take 2 cm<sup>3</sup> of urine to glass tube, acidify with few drops of 2 mol/dm<sup>3</sup> HNO<sub>3</sub> and add few drops of 0.2% AgNO<sub>3</sub> - as positive result white, caseous precipitate of AgCl will be formed.



## - 2. Detection of Phosphate ions

Take 1 cm<sup>3</sup> of urine to glass tube, add 1 cm<sup>3</sup> concentrated HNO<sub>3</sub> and 2 cm<sup>3</sup> ammonium molybdate solution (molibdenian amonu). Heat it carefully over the burner to boil. As positive result yellow precipitate of ammonium phosphoromolybdate will be formed.

## 3. Detection of Calcium ions

Take 2 cm<sup>3</sup> of urine to glass tube, add few drops of 0,2 mol/dm<sup>3</sup> CH<sub>3</sub>COOH and 0,1 mol/dm<sup>3</sup> ammonium oxalate ((NH<sub>4</sub>)<sub>2</sub>C<sub>2</sub>O<sub>4</sub>). As positive result white cloudiness of calcium oxalate (CaC<sub>2</sub>O<sub>4</sub>) will be formed.

## 4. Detection of Sulphate ions

Take 5 cm<sup>3</sup> of urine to glass tube and add few drops of 2 mol/dm<sup>3</sup> HCl and 1 cm<sup>3</sup> of 0,25 mol/dm<sup>3</sup> BaCl<sub>2</sub>. As positive result white cloudiness of BaSO<sub>4</sub> will be formed.

## 3. Properties of pathological urine

### Protocol.

#### 1. Detection of protein

- Take 1 cm<sup>3</sup> of urine and add 2 drops of 10% sulphosalicylic acid. Precipitate or cloudiness will indicate the presence of protein.
- Take 1 cm<sup>3</sup> of urine and add few drops of 10% TriChlorAcetic acid. Cloudiness will indicate the presence of protein.
- Reaction of Heller - take 1 cm<sup>3</sup> of concentrated HNO<sub>3</sub> and gently avoiding mixing of solutions add 1 cm<sup>3</sup> of urine. Thin layer of denatured protein will appear on the border between solutions.

#### 2. Detection of glucose

Take 1 cm<sup>3</sup> of Benedict reagent and add 0,5 cm<sup>3</sup> of urine. Heat it carefully over the burner to boil. As positive result yellow or red precipitate, depending on the content of sugars, will be formed.

#### 3. Detection of ketone bodies - reaction of Legal

Take 1 cm<sup>3</sup> of urine, add 1 cm<sup>3</sup> of 3% nitropruside and 0.5 cm<sup>3</sup> 2 mol/dm<sup>3</sup> NaOH. Red color will be formed.

#### 4. Detection of bilirubin

- Reaction of Gmelin - take 2 cm<sup>3</sup> of urine and add gently 1 cm<sup>3</sup> of concentrated HNO<sub>3</sub>. Green layer on the border between solutions will be formed.
- Reaction of Rosin - take 2 cm<sup>3</sup> of urine and add gently 0,05 % solution of iodine. Grey-green colour caused by the oxidation of bilirubin to biliwerdin will be formed.

