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LITERATURA

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2):116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

Data wydania: 10. 2023.

A-400 BIL TOTAL

Cat. No **7-445** (EN)

INTENDED USE

Diagnostic kit for determination of total bilirubin concentration intended to use in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions

INTRODUCTION

Bilirubin is a yellow pigment – product of heme degradation. For clinical purposes, bilirubin is expressed as two fractions: conjugated and unconjugated. In hepatocytes bilirubin is enzymatically conjugated with glucuronic acid residues. This form is called direct or conjugated. Bilirubin without glucuronic acid modification is bound to albumin and is termed unconjugated or indirect. Indirect bilirubin is calculated as the difference between total and direct bilirubin.

Serum bilirubin measurement is widely used as a screening test for liver functions. Hiperbilirubinemia is usually the result of jaundice (mechanical, hemolytic), Dubin-Johnson syndrome, Gilbert's syndrome, Crigler-Najjar syndrome, bile ducts disease.

METHOD PRINCIPLE

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and vanadate in an acidic solution, total bilirubin (both conjugated – direct, and unconjugated bilirubin) is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the total bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

REAGENTS

Package

1-REAGENT 2 x 31 ml
2-REAGENT 1 x 17 ml

The reagents when stored at 10-25°C are stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

Concentrations in the test

1-REAGENT
citrate buffer (pH 2.8) 90 mmol/l
detergent
2-REAGENT
phosphate buffer (pH 7.0) 4.6 mmol/l
sodium metavanadate 3.0 mmol/l

Warnings and notes

- Do not freeze reagents.
- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not interchange caps.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H319 Causes serious eye irritation.
H410 Very toxic to aquatic life with long lasting effects.

P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P273 Avoid release to the environment.

P391 Collect spillage.

SPECIMEN

Serum free from haemolysis.

Serum should be separated from red blood cells as soon as possible after blood collection.

Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended. It is recommended to follow CLSI procedures regarding specimen collecting and handling.

Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark at 2-8°C, at the most 3 days.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-REAGENT and 2-REAGENT are ready to use.

For reagent blank deionized water is recommended.

Actions required:

When performing assays at analyzer **BS-400**, there is a probability of **cross-contamination** affecting the tests results: CREATININE – BIL TOTAL II GEN. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER.

REFERENCE VALUES⁴

Age	mg/dl	μmol/l
0 - 1 days	< 8	< 137
1 - 2 days	< 12	< 205
3 - 5 days	< 16	< 274
5 days - 60 years	0.3 - 1.2	5 - 21
60 - 90 years	0.2 - 1.1	3 - 19
> 90 years	0.2 - 0.9	3 - 15

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems are recommended the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177). As a 0 calibrator deionized water should be used.

The calibration curve should be prepared every 9 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity:

0.06 mg/dl (1.03 μmol/l) – BS-400

0.07 mg/dl (1.20 μmol/l) – BS-480

Linearity:

up to 69 mg/dl (1180 μmol/l) – BS-400

up to 48.7 mg/dl (832.9 μmol/l) – BS-480

Specificity / Interferences

Haemoglobin up to 0.25 g/dl, ascorbic acid up to 500 mg/l and intralipid up to 250 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	level 1	1.00	0.02	1.73
	level 2	4.20	0.02	0.49
BS-480 (n = 10)	level 1	1.17	0.01	0.49
	level 2	4.33	0.01	0.13
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 19)	level 1	1.03	0.08	7.26
	level 2	4.20	0.14	3.28
BS-480 (n = 10)	level 1	1.20	0.01	1.20
	level 2	4.46	0.05	1.21

Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 41 samples gave following results:

$y = 0.8872x + 0.0143$ mg/dl;

$R = 0.999$ (R – correlation coefficient)

A comparison between total bilirubin values determined at **BS-480** (y) and at **Advia 1650** (x) using 82 samples gave following results:

$y = 0.9187x + 0.0385$ U/l;

$R = 0.999$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Tokuda K, Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2):116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2);116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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• Calibration		Judgment Criteria			
Calibration					
Rule	<input type="text" value="Multi-point Linear"/>	Sensitivity	<input type="text"/>	Blank Abs.	<input type="text"/>
Replicate	<input type="text" value="3"/>	Factor Diff.	<input type="text"/>	Error Limit	<input type="text"/>
K	<input type="text"/>	SD	<input type="text"/>	Corr. Coeff.	<input type="text"/>
• QC				Auto QC	
Rules					
Westgard Multi-rule		Cum. Sum Check		Interval <input type="text"/>	
<input type="checkbox"/> 1-2S	<input type="checkbox"/> R-4S	<input type="checkbox"/> 1.0 - 2.7			
<input type="checkbox"/> 1-3S	<input type="checkbox"/> 4-1S	<input checked="" type="checkbox"/> 1.0 - 3.0			
<input type="checkbox"/> 2-2S	<input type="checkbox"/> 10-X	<input type="checkbox"/> 0.5 - 5.1			

• BS-480

Chem	<input type="text" value="BIL T"/>	No.	<input type="text" value="005"/>	Sample Type	<input type="text" value="SERUM"/>
Chemistry	<input type="text" value="Bilirubin Total"/>	Print name	<input type="text" value="BIL T"/>		
Reaction Type	<input type="text" value="Endpoint"/>	Reaction Direction	<input type="text" value="Decrease"/>		
Pri Wave	<input type="text" value="450"/>	Sec Wave	<input type="text" value="546"/>		
Unit	<input type="text" value="mg/dl"/>	Decimal	<input type="text" value="0.01"/>		
Blank Time	<input type="text" value="46"/> <input type="text" value="47"/>	Reaction Time	<input type="text" value="81"/> <input type="text" value="82"/>		
Standard	Sample Vol <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/>	Aspirated	<input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/> <input type="text" value="20"/>	Diluent	<input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/> <input type="text" value="180"/>
Decreased	<input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/>	Reagent Vol	<input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/>	Diluent	<input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/>
Increased	<input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/>	R1	<input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/> <input type="text" value="240"/>	R2	<input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/>
	<input type="checkbox"/> Sample Blank <input type="checkbox"/> Auto Retun	R3	<input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/>	R4	<input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/> <input type="text" value="60"/>
Linearity Range (Standard)	<input type="text" value="0.07"/> <input type="text" value="48.7"/>	Linearity Limit	<input type="text"/>		
Linearity Range (Decreased)	<input type="text"/>	Substrate Depletion	<input type="text"/>		
Linearity Range (Increased)	<input type="text"/>	Mixed Blank Abs	<input type="text" value="-33000"/> <input type="text" value="33000"/>		
R1 Blank Abs	<input type="text" value="-33000"/> <input type="text" value="33000"/>	Uncapping Time	<input type="text" value="84"/> Day(s)		
Blank Response	<input type="text" value="-33000"/> <input type="text" value="33000"/>	Reagent Alarm Limit	<input type="text"/>		
Twin Chemistry	<input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check	• Antigen Addition			
Q1 <input type="text" value="0"/>	Q2 <input type="text" value="0"/>	Q3 <input type="text" value="0"/>	Q4 <input type="text" value="0"/>		
PC <input type="text" value="0"/>	ABS <input type="text" value="0"/>				

Calibration Settings		Auto Calibration	
Math Model	<input type="text" value="Multi-point Linear"/>	<input type="checkbox"/> Bottle Changed	
Factor	<input type="text"/>	<input type="checkbox"/> Lot Changed	
Replicates	<input type="text" value="3"/>	<input type="checkbox"/> Cal Time	
Acceptance Limits			
Cal Time	<input type="text" value="2016"/> Hour	SD	<input type="text"/>
Slope Diff	<input type="text"/>	Repeatability	<input type="text"/>
Sensitivity	<input type="text"/>		
Deter Coeff	<input type="text"/>		

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LITERATURA

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24. 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Praventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembińska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

Data wydania: 06. 2023.

A-400 CHOL

Cat. No **7-404**

(EN)

Warning

H319 Causes serious eye irritation.
P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis. Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for cholesterol measurement. Plasma cholesterol values have been reported to be 3% to 5% lower than serum cholesterol values. Serum should be separated from red blood cells as soon as possible after blood collection. Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

I-REAGENT is ready to use.
For reagent blank deionized water is recommended.

Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: CHOL – LIPASE II GEN, CHOL – LIPASE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER.

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QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2 (Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analyzer system BS-400 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and/or LEVEL 2 (Cat. No 5-175; 5-177). **These two calibrators or calibrator and deionised water** should be used for calibration.

For the calibration of automatic analyzer system BS-480 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity:

20 mg/dl (0.518 mmol/l) – BS-400
12 mg/dl (0.311 mmol/l) – BS-480

Linearity:

up to 750 mg/dl (19.4 mmol/l) – BS-400
up to 800 mg/dl (20.7 mmol/l) – BS-480

Specificity / Interferences

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	level 1	156.11	0.62	0.40
	level 2	293.39	1.88	0.64
BS-480 (n = 10)	level 1	93.45	0.31	0.33
	level 2	244.81	0.81	0.33
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 14)	level 1	158.91	1.41	0.89
	level 2	298.05	2.24	0.75
BS-480 (n = 10)	level 1	92.55	0.67	0.72
	level 2	239.64	1.48	0.62

Method comparison

A comparison between cholesterol values determined at **BS-400** (y) and at **Olympus AU640** (x) using 40 samples gave following results:

$y = 0.9112x + 8.1417$ mg/dl;
 $R = 0.996$ (R – correlation coefficient)

A comparison between cholesterol values determined at **BS-480** (y) and at **Cobas Integra 400 Plus** (x) using 40 samples gave following results:

$y = 0.9826x + 3.3927$ mg/dl;
 $R = 0.993$ (R – correlation coefficient)

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LITERATURE

- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29. 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18. 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24. 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Präventivmed. 10. 25 (1975).
- Richmond W.: Clin. Chem. 19. 1350 (1973).
- Roeschlau P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12. 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia. PA: WB Saunders. 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6. 24 (1969).
- Dembińska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts) Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148). 36-69. (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed.. Laboratory Tests Handbook. 4th ed.. Hudson. Lexi-Comp. 143. (1996).

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- Siedel J., Hägele E.O., Zigenhorn J., Wahlefeld A.W.: Clin. Chem. 29, 6 (1983).
- Tel R.M., Berends G.T.: J Clin.Chem. Clin. Biochem. 18, 10 (1980).
- Rautela G.S., Liedtke R.J.: Clin. Chem. 24, 1 (1978)
- Schettler G., Nussel E.: Arbeitsmed. Sozialmed. Praventivmed. 10, 25 (1975).
- Richmond W.: Clin. Chem. 19, 1350 (1973).
- Roeschlaup P., Bernt E., Gruber W.: J. Clin. Chem. Biochem. 12, 403 (1974).
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA: WB Saunders, 130 (1995).
- Trinder P.: Ann. Clin. Biochem. 6, 24 (1969).
- Dembińska-Kieć A., Naskalski J.W., Solnica B.: Diagnostyka laboratoryjna z elementami biochemii klinicznej. Wyd. IV, str. 670, 875 (2017).
- Cybulska B, Szostak WB, Filipiak KJ, et al. Polish Forum for Prevention Guidelines on Dyslipidaemia: update 2016. Kardiol Pol. 2017; 75(2): 187–190, doi: 10.5603/KP.2017.0031.
- Piepoli MF, Hoes AW, Agewall S, et al. 2016 European Guidelines on cardiovascular disease prevention in clinical practice: The Sixth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of 10 societies and by invited experts)Developed with the special contribution of the European Association for Cardiovascular Prevention & Rehabilitation (EACPR). Eur Heart J. 2016; 37(29): 2315–2381, doi: 10.1093/eurheartj/ehw106, indexed in Pubmed: 27222591.
- NCEP Expert Panel. Arch Inter Med (148), 36-69, (1988).
- Jacobso D.S., DeMott W.R., Grady H.J., et. al., ed., Laboratory Tests Handbook, 4th ed., Hudson, Lexi-Comp, 143, (1996).

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Calibration

Calibration	Judgment Criteria		
Rule	Two-point Linear	Sensitivity	
Replicate	3	Factor Diff.	
K		SD	
QC		Blank Abs.	
Rules		Error Limit	
Westgard Multi-rule		Corr. Coeff.	
v 1-2S	v R-4S	Cum. Sum Check	1.0 - 2.7
v 1-3S	v 4-1S		1.0 - 3.0
v 2-2S	v 10-X		0.5 - 5.1
		Auto QC	Interval

BS-480

Chem	CHOL	No.	002	Sample Type	SERUM
Chemistry	CHOLESTEROL			Print name	CHOL
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	505	Sec Wave	800		
Unit	mg/dL	Decimal	1		
Blank Time	11 12	Reaction Time	44 46		
Standard	Sample Vol 2.5 μL Aspirated 20 μL Diluent 180 μL	Reagent Vol	200 μL	Diluent	
Decreased	2.5 μL 20 μL 180 μL	R1			
Increased		R2			
		R3			
		R4			
Linearity Range	12 800	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	77 Dav(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension			
	Prozone Check o Rate Check	Antigen Addition			
Q1 0 Q2 0 Q3 0 Q4 0					
PC 0 ABS 0					

Calibration Settings

Math Model	Multi-point Linear	Bottle Changed	
Factor		Lot Changed	
Replicates	3	Cal Time	

Acceptance Limits

Cal Time	1848 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			